




Source of information: leaflet (note: the instructions might differ according to the brand of the acquired inhaler)

Unpack your inhaler

The inhaler comes with a laminate pouch to keep the powder dry. Do not open the pouch until you are ready to start using your inhaler.


Using your inhaler

1



Open the protective cover, if you are using it. Remove the dust cap.

2



Shake it vigorously up and down three to five times. Don't press the top while you shake it.

3



Press the top of the inhaler once between your finger and thumb until you hear a click, and then release it.

4



Breathe out as far as is comfortable. Do not breathe in again yet.

5



Place the mouthpiece in your mouth between your teeth. Close your lips around it. Do not bite.

6



Take a strong and deep breath through your mouth. Remove the Easyhaler from your mouth.

7



Hold your breath for at least 5 seconds.

8



Breathe out gently.

9



Put the dust cap back on the mouthpiece. Close the protective cover, if you are using it.

