

## my-inhaler.com


forspiro

Source of information: leaflet (note: the instructions might differ according to the brand of the acquired inhaler)

## Before using the inhaler

The transparent side chamber door should be opened. The foil strip should be removed from the side chamber by carefully tearing away the full length of strip against the "teeth" of the side chamber. The strip should not be pulled or tugged. The side chamber door should be closed and the used strip should be disposed of. Note: As the inhaler is used the side chamber will gradually fill up with used strip. The foil strips with black bars don't contain medication. Eventually the numbered sections of the strip will appear in the side chamber. There should never be more than 2 sections of foil strip in the side chamber as they may cause the inhaler to jam. The strip should be torn away carefully, and disposed of safely. Also remember: The white lever should only be operated when the patient is ready to inhale a dose of the medication. If the patient plays with the white lever he/she will waste doses.

## Using your inhaler





my-inhaler.com

forspiro



The inhaler should be removed from the mouth and the breath should be held for 5-10 seconds or as long as is possible without causing discomfort.



Afterwards, the patient should breathe out slowly, but not into the inhaler.



The protective cap should be closed over the mouthpiece.

